REAL 6100 PLUS ADULT

Order form - your support for assessment, quotation and order

Use *Adobe Acrobat Reader* to complete this form digitally. Download the program **here**.

☐ Quotation ☐ Order Number of chairs _____

YOUR CONTACT INFORMATION

Name _____ Email _____ Phone _____

1 MODEL

Power module

☐ LiNX

Color of seat, backrest and other upholstery (except armrests):

Plush
☐ Grey ☐ REAL-blue
☐ Bordeaux ☐ Black

Artificial leather
☐ Red ☐ Black
☐ Grey ☐ Blue
☐ Brown

Hygiene
☐ Grey

Patterned
☐ Blue

2 SEAT

SEAT UPHOLSTERY

Width x depth, in cm	40x39	40x43	40x48	44x43	44x48	44x53	48x48	48x53
ErgoMedic (EM)	☐	☐	☐	☐	☐	☐	☐	☐
EM SlowRecovery	☐	☐	☐	☐	☐	☐	☐	☐
EM PLUS	☐	☐	☐	☐	☐	☐	☐	☐
EM PLUS SlowRecovery	☐	☐	☐	☐	☐	☐	☐	☐
EM wooden seat plate	☐	☐	☐	☐	☐	☐	☐	☐
EM wooden seat plate with Anti-slip	☐	☐	☐	☐	☐	☐	☐	☐

Matching tips!
(measurements refer to the width of the upholstery, in cm)

Seat	Backrest
40	33
44	38-42
48	38-46

SEAT TILT seat and backrest are tilted together [1]

Fixed interval
☐ No tilt
☐ 15° fw, 8° bw
☐ 15° fw, 15° bw

Allangle (= configurable)
☐ 15° fw, 8° bw
☐ 8° fw, 15° bw
☐ 0° fw, 23° bw [2]

Seat angle control
☐ Handle control under seat
☐ Gas spring at armrest
☐ Electrical

Armrests follow tilt [3]
☐ Yes
☐ No

[1] *See section 7*
[2] *Requires rear base extenders*
[3] *Only for seat tilt "Allangle" and "15° fw, 15° bw"*

3 BACKREST

BACKREST UPHOLSTERY	Low				High				
Width x height, in cm	33x26	35x20	38x26	40x20	33x36	35x43	40x34	42x43	46x43
ErgoMedic (EM) [4]	☐		☐		☐		☐	☐	☐
EM SlowRecovery	☐		☐		☐		☐	☐	
EM Support						☐		☐	☐
EM PLUS		☐		☐		☐		☐	☐
EM PLUS SlowRecovery		☐		☐		☐		☐	☐
EM PLUS side wedge						☐		☐	☐
EM PLUS Comfort (= lower back support)						☐		☐	☐
EM PLUS side wedge Comfort						☐		☐	☐
EM PLUS reinforced side wedge						☐		☐	☐
EM PLUS reinforced side wedge Comfort						☐		☐	☐
EM PLUS without bulb								☐	☐
EM PLUS without bulb w side wedge								☐	☐

[4] *46x43 only available in grey plush*

(more backrests on next page)

BACKREST UPHOLSTERY (CONTINUATION)	Low		High	
	39x30	45x30	39x40	45x40
Slingback (backrest with adjustable straps)	☐	☐	☐	☐

BACKREST MECHANISM

Model

- ☐ Standard
- ☐ Medic
- ☐ Double Medic (Allangle required)
- ☐ Comfort
- ☐ Electric
- ☐ Slingback

Control position (only for comfort)

- ☐ Armrest right front
- ☐ Armrest left front
- ☐ Armrest right back
- ☐ Armrest left back

4 ARMREST

ARMREST PLATE black molded / artificial leather

Length in cm	25	30	35	40	42	45
Mobile	☐	☐	☐	☐		
ErgoMedic	☐	☐	☐	☐		
Extra Soft	☐	☐	☐			
ErgoMedic PLUS		L ☐☐ R			L ☐☐ R	
Hemiplegia						L ☐☐ R

ARMREST MECHANISM

Model

- ☐ PLUS (non foldable)
- ☐ Foldable backwards

5 CHASSIS

ACTUATOR FOR SEAT HEIGHT

Electric seat height adjustment, in cm [5]

- ☐ Low 38-58
- ☐ High 46-74

User weight

- ☐ 0-135 kg [6]

[5] *Seat height measured between floor and bottom of seat (without seat tilt)*
[6] *For user weight exceeding 135 kg, contact Mercado customer service*

CHASSIS (6100 PLUS)

Model

- ☐ 6-wheels with LiNX electrical system

Rear leg extenders

- ☐ Yes ☐ No

6 OPTIONAL AND ALTERNATIVE EQUIPMENT

JOYSTICK

Model

- ☐ LiNX REM211
- ☐ LiNX REM400
- ☐ LiNX REM500

Position

- ☐ Right armrest
- ☐ Left armrest
- ☐ On table
- ☐ Attendant remote control

Attachment

- ☐ ParaLoc
- ☐ Operating arm long (115 mm)
- ☐ Operating arm short (95 mm)

Primary controls for seating functions

- ☐ With joystick
- ☐ With 10 button keypad LiNX

Special joystick

- ☐ Mini-joystick [7]
- ☐ Attendant control
- ☐ Head control [7]
- ☐ Heavy duty joystick [7]

Alternative joystick shapes

- ☐ Half-moon Ø50 mm
- ☐ Ball Ø30 mm
- ☐ Ball Ø50 mm
- ☐ Stick Ø30 x h80 mm
- ☐ Soft yellow ball Ø30 mm

[7] *Not compatible with REM211*

Mercado Medic AB
Tryffelslingan 14, S-181 57 Lidingö
Sweden

Tel 08 555 143 00

Email kundtjanst@mercado.se
Web www.mercado.se

LEGREST Adjustable interval, in cm

Legrest Medic
- ☐ M 23-32
- ☐ L 28-41
- ☐ XL 38-50

Legrest Medic divided
- ☐ M 23-32
- ☐ L 28-41
- ☐ XL 38-50

Legrest center mounted
- ☐ 20-31 (W29xD24)
- ☐ 30-38 (W29xD24)
- ☐ 37-53 (W29xD24)
- ☐ 20-31 Mini (W25xD18)
- ☐ 30-38 Mini (W25xD18)
- ☐ 37-53 Mini (W25xD18)
- ☐ 30-38 Slim (W32xD10)
- ☐ 37-53 Slim (W32xD10)

- ☐ Rim on back and sides of footplate
- ☐ Electrical angle legrest (0-90°)
- ☐ Electric length compensating legrest (0-90°)

Adapter for Cross legrest
- ☐ Right
- ☐ Left

- ☐ Including Cross legrest
- ☐ 23-31 cm
- ☐ 33-51 cm

EXTRA SUPPORT Width x height, in cm

Trunk support [8]
- ☐ Straight
- ☐ Curved
- ☐ Curved and adjustable

Side support
- ☐ 28x8
- ☐ 24x14
- ☐ 30x16

Hip belt (length)
- ☐ M 63-83
- ☐ L 83-103

Attachment for belt/harness
- ☐ Attachment hip belt [9]
- ☐ Attachment H-harness

Neckrest [8]
- ☐ Flat 29x16
- ☐ Rounded 31x18
- ☐ ErgoMedic PLUS 30x23
- ☐ ErgoMedic PLUS 30x8

Neckrest mechanism [8]
- ☐ Height adjustable
- ☐ Height and sideways adjustable

Knee support
- ☐ 16x8

Calf support
- ☐ Medic divided
- ☐ Center mounted

- ☐ 10x10
- ☐ 15x13

Divider
- ☐ Leg divider
- ☐ Knee divider

- ☐ 11 (width)

[8] *Only for high backrest*
[9] *Not needed if the selected option for seat tilt is "Allangle"*

BATTERY CHARGER
- ☐ Yes
- ☐ No [10]

UPTIME METER
- ☐ Yes
- ☐ No

COVER
- ☐ Backrest, plush
- ☐ Seat, plush
- ☐ Backrest, water resistant
- ☐ Seat, water resistant

IMPACT PROTECTION
- ☐ Impact protection for front wheels

TABLE/TRAY
- ☐ Table 37x45 cm foldable
- ☐ Table 60x45 cm foldable
- ☐ Tray 45x30 cm

[10] *CE marking not valid if product is used with charger other than chargers tested and approved by Mercado Medic*

7 PRESCRIPTION RECOMMENDATIONS

The backrest mechanism Double Medic should be prescribed when:

... there is a risk for involuntary movements which may cause great stress on the backrest mechanism over a long period of time.
... the user is over 200 cm long and has a backward tilted pelvis, in combination with a seat unit that can tilt more than 8° backward, and neck support.
... the seat unit can be tilted more than 15° backward and the user weighs more than 100 kg.
... the reinforced seat frame Allangle +15°/-23° is selected due to high loads.
... a user with spasms or repeated rocking movements is going to use the chair.
... backrest upholstery wider than 46 cm is mounted on the chair in combination with 15° tilt backward.

COMMENTS